


Product Datasheet

C7orf16 (PPP1R17) antibody (orb1315105)

Description	Mouse monoclonal antibody to C7orf16 (PPP1R17)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human PPP1R17(NP_006649) produced in E.coli.
Target	PPP1R17
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	17.7 kDa
Uniprot ID	O96001
NCBI	NM_006658
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...