


Product Datasheet

EBLN2 antibody (orb1315143)

Description	Mouse monoclonal antibody to EBLN2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human EBLN2(NP_060499) produced in E.coli.
Target	EBLN2
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	30.3 kDa
Uniprot ID	Q6P217
NCBI	NM_018029
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...