


Product Datasheet

SNX10 antibody (orb1315164)

Description	Mouse monoclonal antibody to SNX10
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human SNX10(NP_037454) produced in E.coli.
Target	SNX10
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	23.4 kDa
Uniprot ID	Q9Y5X0
NCBI	NM_013322
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...