


Product Datasheet

C19orf80 (ANGPTL8) antibody (orb1315170)

Description	Mouse monoclonal antibody to C19orf80 (ANGPTL8)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Human recombinant protein fragment corresponding to amino acids 22-198 of human C19ORF80(NP_061157) produced in E.coli.
Target	ANGPTL8
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	21.9 kDa
Uniprot ID	Q6UXH0
NCBI	NM_018687
Dilution Range	IHC 1:150
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...