


Product Datasheet

c Kit (KIT) antibody (orb1315251)

Description	Mouse monoclonal antibody to c Kit (KIT)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 27-308 of human KIT(NP_000213) produced in E.coli.
Target	KIT
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	109.7 kDa
Uniprot ID	P10721
NCBI	NM_000222
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...