


Product Datasheet

Her2 (ERBB2) antibody (orb1315269)

Description	Mouse monoclonal antibody to Her2 (ERBB2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	HER2 extracellular domain (aa 23-652) expressed in HEK293 cells with a C-terminal DDK tag
Target	ERBB2
Note	For research use only
Isotype	IgG2b, kappa
Clonality	Monoclonal
Uniprot ID	P04626
NCBI	NM_004448
Dilution Range	FC: 1:50
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of live MCF7 cells...