


Product Datasheet

HPCAL4 antibody (orb1315278)

Description	Mouse monoclonal antibody to HPCAL4
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 2-191 of human HPCAL4(NP_057341) produced in E.coli.
Target	HPCAL4
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	22 kDa
Uniprot ID	Q9UM19
NCBI	NM_016257
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...