


## Product Datasheet

### C19orf80 (ANGPTL8) antibody (orb1315283)

|                            |                                                                                                                         |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to C19orf80 (ANGPTL8)                                                                         |
| <b>Species/Host</b>        | Mouse                                                                                                                   |
| <b>Reactivity</b>          | Human                                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                                            |
| <b>Tested Applications</b> | IHC, WB                                                                                                                 |
| <b>Immunogen</b>           | Human recombinant protein fragment corresponding to amino acids 22-198 of human C19ORF80(NP_061157) produced in E.coli. |
| <b>Target</b>              | ANGPTL8                                                                                                                 |
| <b>Note</b>                | For research use only                                                                                                   |
| <b>Isotype</b>             | IgG1                                                                                                                    |
| <b>Clonality</b>           | Monoclonal                                                                                                              |
| <b>MW</b>                  | 21.9 kDa                                                                                                                |
| <b>Uniprot ID</b>          | <a href="#">Q6UXH0</a>                                                                                                  |
| <b>NCBI</b>                | <a href="#">NM_018687</a>                                                                                               |
| <b>Dilution Range</b>      | WB 1:2000, IHC 1:150                                                                                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                         |



Western blot analysis of HEK293T cells w...