


Product Datasheet

GAS8 antibody (orb1315313)

Description	Mouse monoclonal antibody to GAS8
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human GAS8 (NP_001472) produced in E.coli.
Target	GAS8
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	56.2 kDa
Uniprot ID	O95995
NCBI	NM_001481
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...