


Product Datasheet

ZNF449 antibody (orb1315331)

Description	Mouse monoclonal antibody to ZNF449
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human ZNF449 (NP_689908) produced in E.coli.
Target	ZNF449
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	59.8 kDa
Uniprot ID	Q6P9G9
NCBI	NM_152695
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...