


Product Datasheet

NAP1L3 antibody (orb1315369)

Description	Mouse monoclonal antibody to NAP1L3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-315 of human NAP1L3(NP_004529) produced in E.coli.
Target	NAP1L3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	57.4 kDa
Uniprot ID	Q99457
NCBI	NM_004538
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...