


Product Datasheet

CHMP5 antibody (orb1315391)

Description	Mouse monoclonal antibody to CHMP5
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human CHMP5 (NP_057494) produced in E.coli.
Target	CHMP5
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	24.4 kDa
Uniprot ID	Q9NZZ3
NCBI	NM_016410
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...