


Product Datasheet

EXOSC1 antibody (orb1315401)

Description	Mouse monoclonal antibody to EXOSC1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human EXOSC1 (NP_057130) produced in E.coli.
Target	EXOSC1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	21.3 kDa
Uniprot ID	Q9Y3B2
NCBI	NM_016046
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...