


Product Datasheet

IHO1 antibody (orb1315425)

Description	Mouse monoclonal antibody to IHO1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CCDC36(NP_835467) produced in E.coli.
Target	IHO1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	64.9 kDa
Uniprot ID	Q8IYA8
NCBI	NM_178173
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...