


Product Datasheet

C2orf80 antibody (orb1315437)

Description	Mouse monoclonal antibody to C2orf80
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human C2orf80 (NP_001092804) produced in E.coli.
Target	C2orf80
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q0P641
NCBI	NM_001099334
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...