


Product Datasheet

TTF1 (NKX2-1) antibody (orb1315444)

Description	Mouse monoclonal antibody to TTF1 (NKX2-1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human NKX2-1(NP_003308) produced in E.coli.
Target	NKX2-1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	38.4 kDa
Uniprot ID	P43699
NCBI	NM_003317
Dilution Range	IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...