


Product Datasheet

SPATA22 antibody (orb1315455)

Description	Mouse monoclonal antibody to SPATA22
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human SPATA22(NP_115987) produced in E.coli.
Target	SPATA22
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	41.1 kDa
Uniprot ID	Q8NHS9
NCBI	NM_032598
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...