

Product Datasheet

MAP3K8 antibody (orb1315482)

Description	Mouse monoclonal antibody to MAP3K8
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 277-461 of human MAP3K8(NP_005195) produced in E.coli.
Target	MAP3K8
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	52.7 kDa
Uniprot ID	P41279
NCBI	NM_005204
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...