


## Product Datasheet

### IHO1 antibody (orb1315484)

|                            |                                                                                      |
|----------------------------|--------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to IHO1                                                    |
| <b>Species/Host</b>        | Mouse                                                                                |
| <b>Reactivity</b>          | Human                                                                                |
| <b>Conjugation</b>         | Unconjugated                                                                         |
| <b>Tested Applications</b> | IHC, WB                                                                              |
| <b>Immunogen</b>           | Full length human recombinant protein of human CCDC36(NP_835467) produced in E.coli. |
| <b>Target</b>              | IHO1                                                                                 |
| <b>Note</b>                | For research use only                                                                |
| <b>Isotype</b>             | IgG1                                                                                 |
| <b>Clonality</b>           | Monoclonal                                                                           |
| <b>MW</b>                  | 64.9 kDa                                                                             |
| <b>Uniprot ID</b>          | <a href="#">Q8IYA8</a>                                                               |
| <b>NCBI</b>                | <a href="#">NM_178173</a>                                                            |
| <b>Dilution Range</b>      | WB 1:2000, IHC 1:150                                                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                      |



Western blot analysis of HEK293T cells w...