


Product Datasheet

Cytokeratin 16 (KRT16) antibody (orb1315485)

Description	Mouse monoclonal antibody to Cytokeratin 16 (KRT16)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human KRT16 (NP_005548) produced in E.coli.
Target	KRT16
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	51.1 kDa
Uniprot ID	P08779
NCBI	NM_005557
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...