


Product Datasheet

RAB11FIP4 antibody (orb1315490)

Description	Mouse monoclonal antibody to RAB11FIP4
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 91-177 of human RAB11FIP4(NP_116321) produced in E.coli.
Target	RAB11FIP4
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	71.7 kDa
Uniprot ID	Q86YS3
NCBI	NM_032932
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...