


Product Datasheet

FMO3 antibody (orb1315536)

Description	Mouse monoclonal antibody to FMO3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-291 of human FMO3(NP_008825) produced in E.coli.
Target	FMO3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P31513
NCBI	NM_006894
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...