


Product Datasheet

Telethonin (TCAP) antibody (orb1315555)

Description	Mouse monoclonal antibody to Telethonin (TCAP)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human TCAP (NP_003664) produced in E.coli.
Target	TCAP
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	18.9 kDa
Uniprot ID	O15273
NCBI	NM_003673
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...