


Product Datasheet

C4orf36 antibody (orb1315581)

Description	Mouse monoclonal antibody to C4orf36
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human C4orf36(NP_653246) produced in E.coli.
Target	C4orf36
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	13.2 kDa
Uniprot ID	Q96KX1
NCBI	NM_144645
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...