


Product Datasheet

TBC1D28 antibody (orb1315598)

Description	Mouse monoclonal antibody to TBC1D28
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human TBC1D28 (NP_001034486) produced in E.coli.
Target	TBC1D28
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	23.9 kDa
Uniprot ID	Q2M2D7
NCBI	NM_001039397
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...