

Product Datasheet

ZFP200 (ZNF200) antibody (orb1315607)

Description	Mouse monoclonal antibody to ZFP200 (ZNF200)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human ZNF200(NP_932353) produced in E.coli.
Target	ZNF200
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	45.2 kDa
Uniprot ID	P98182
NCBI	NM_198087
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.

