


Product Datasheet

PSMG2 antibody (orb1315635)

Description	Mouse monoclonal antibody to PSMG2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human PSMG2 (NP_064617) produced in E.coli.
Target	PSMG2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	29.2 kDa
Uniprot ID	Q969U7
NCBI	NM_020232
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...