


Product Datasheet

SNRPB2 antibody (orb1315653)

Description	Mouse monoclonal antibody to SNRPB2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human SNRPB2 (NP_937863) produced in E.coli.
Target	SNRPB2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	25.3 kDa
Uniprot ID	P08579
NCBI	NM_198220
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...