


Product Datasheet

IL12A antibody (orb1315659)

Description	Mouse monoclonal antibody to IL12A
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 57-253 of human IL12A(NP_000873) produced in E.coli.
Target	IL12A
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P29459
NCBI	NM_000882
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...