


Product Datasheet

SPT3 (SUPT3H) antibody (orb1315665)

Description	Mouse monoclonal antibody to SPT3 (SUPT3H)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human SUPT3H(NP_852001) produced in E.coli.
Target	SUPT3H
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	O75486
NCBI	NM_181356
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...