


Product Datasheet

SDHB antibody (orb1315668)

Description	Mouse monoclonal antibody to SDHB
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human SDHB(NP_002991) produced in E.coli.
Target	SDHB
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	28.8 kDa
Uniprot ID	P21912
NCBI	NM_003000
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...