

Product Datasheet

MS4A12 antibody (orb1315676)

Description	Mouse monoclonal antibody to MS4A12
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human MS4A12 (NP_060186) produced in HEK293T cell.
Target	MS4A12
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	27.9 kDa
Uniprot ID	Q9NXJ0
NCBI	NM_017716
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.

