

Product Datasheet

Cardiac Troponin T (TNNT2) antibody (orb1315716)

Description Mouse monoclonal antibody to Cardiac Troponin T (TNNT2)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC, WB

Immunogen Native troponin complex

Target TNNT2

Note For research use only

Isotype IgG1


Clonality Monoclonal

Uniprot ID [P45379](#)

NCBI [NM_000364](#)

Dilution Range WB 1:2000, IHC 1:150

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...