


Product Datasheet

SP110 antibody (orb1315743)

Description	Mouse monoclonal antibody to SP110
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 232-521 of human SP110(NP_004501) produced in E.coli.
Target	SP110
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	61.7 kDa
Uniprot ID	Q9HB58
NCBI	NM_004510
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...