


Product Datasheet

IL3RB (CSF2RB) antibody (orb1315745)

Description	Mouse monoclonal antibody to IL3RB (CSF2RB)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human CSF2RB (NP_000386) produced in HEK293T cell.
Target	CSF2RB
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	95.7 kDa
Uniprot ID	P32927
NCBI	NM_000395
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...