


Product Datasheet

SOX21 antibody (orb1315798)

Description	Mouse monoclonal antibody to SOX21
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-140 of human SOX21(NP_009015) produced in E.coli.
Target	SOX21
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q9Y651
NCBI	NM_007084
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...