


Product Datasheet

GPNMB antibody (orb1315854)

Description	Mouse monoclonal antibody to GPNMB
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human GPNMB (NP_001005340) produced in HEK293T cell.
Target	GPNMB
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	61.5 kDa
Uniprot ID	Q14956
NCBI	NM_001005340
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...