


Product Datasheet

GNMB antibody (orb1315860)

Description	Mouse monoclonal antibody to GNMB
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human GNMB (NP_001005340) produced in HEK293T cell.
Target	GNMB
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	Q14956
NCBI	NM_001005340
Dilution Range	WB 1:200~2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...