


Product Datasheet

LIME (LIME1) antibody (orb1315906)

Description	Mouse monoclonal antibody to LIME (LIME1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 30-295 of human LIME1(NP_060276) produced in E.coli.
Target	LIME1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	31.1 kDa
Uniprot ID	Q9H400
NCBI	NM_017806
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...