


Product Datasheet

IKBKE antibody (orb1315930)

Description	Mouse monoclonal antibody to IKBKE
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 474-716 of human IKBKE(NP_054721) produced in E.coli.
Target	IKBKE
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	80.3 kDa
Uniprot ID	Q14164
NCBI	NM_014002
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...