


Product Datasheet

NUP43 antibody (orb1315962)

Description	Mouse monoclonal antibody to NUP43
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-253 of human NUP43(NP_942590) produced in E.coli.
Target	NUP43
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	42 kDa
Uniprot ID	Q8NFH3
NCBI	NM_198887
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...