


Product Datasheet

LOXL2 antibody (orb1316031)

Description	Mouse monoclonal antibody to LOXL2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Human recombinant protein fragment corresponding to amino acids 152-450 of human LOXL2(NP_002309) produced in E.coli.
Target	LOXL2
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	84 kDa
Uniprot ID	Q9Y4K0
NCBI	NM_002318
Dilution Range	IHC 1:150
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...