


Product Datasheet

CD13 (ANPEP) antibody (orb1316049)

Description	Mouse monoclonal antibody to CD13 (ANPEP)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 668-967 of human ANPEP (NP_001141) produced in E.coli.
Target	ANPEP
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P15144
NCBI	NM_001150
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...