


Product Datasheet

SFRP2 antibody (orb1316053)

Description	Mouse monoclonal antibody to SFRP2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 25-295 of human SFRP2(NP_003004) produced in E.coli.
Target	SFRP2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	31.4 kDa
Uniprot ID	Q96HF1
NCBI	NM_003013
Dilution Range	IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...