


Product Datasheet

PON3 antibody (orb1316070)

Description	Mouse monoclonal antibody to PON3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 30-319 of human PON3(NP_000931) produced in E.coli.
Target	PON3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	39.4 kDa
Uniprot ID	Q15166
NCBI	NM_000940
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...