


Product Datasheet

MRPL15 antibody (orb1316078)

Description	Mouse monoclonal antibody to MRPL15
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 28-296 of human MRPL15(NP_054894) produced in E.coli.
Target	MRPL15
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	33.2 kDa
Uniprot ID	Q9P015
NCBI	NM_014175
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...