


Product Datasheet

PISD antibody (orb1316085)

Description	Mouse monoclonal antibody to PISD
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 80-340 of human PISD(NP_055153) produced in E.coli.
Target	PISD
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q9UG56
NCBI	NM_014338
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...