

Product Datasheet

DR5 (TNFRSF10B) antibody (orb1316112)

Description Mouse monoclonal antibody to DR5 (TNFRSF10B)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications IHC

Immunogen Human recombinant protein fragment corresponding to amino acids 203-411 of human TNFRSF10B(NP_671716) produced in E.coli.

Target TNFRSF10B

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 39.4 kDa

Uniprot ID [O14763](#)

NCBI [NM_147187](#)

Dilution Range IHC 1:150

Expiration Date 12 months from date of receipt.



Immunohistochemical staining of paraffin...