


Product Datasheet

NR2C2 antibody (orb1316135)

Description	Mouse monoclonal antibody to NR2C2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-241 of human NR2C2(NP_003289) produced in E.coli.
Target	NR2C2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	67.2 kDa
Uniprot ID	P49116
NCBI	NM_003298
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...