


Product Datasheet

HGF antibody (orb1316145)

Description	Mouse monoclonal antibody to HGF
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 200-494 of human HGF(NP_000592) produced in E.coli.
Target	HGF
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	79.6 kDa
Uniprot ID	P14210
NCBI	NM_000601
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...